


Corrigendum for e-tender

HPLC + 5 Years AMC with UV Detector and Single point control software for Data acquisition and control

S.no.	Component	Specifications in e-tender	Changes in specifications
1	Auto-sampler with sample cooling	Number of Sample Vials: Should accommodate 110 vials or more of 2 ML. Should have advance facility of Auto additions, Priority samples and Auto dilutions.	Number of Sample Vials: Should accommodate 110 vials or more of 1.5 / 2 mL. Auto additions, priority samples, and auto dilutions are desirable features.
2	Column Compartment	Column compartment should have Force Air Circulation technology.	Column compartment should have Force Air Circulation/Block heating and electronic cooling technology.
3	Chromatographic Software	Software to control, acquire and process the data with Oracle data base for high level of data integrity. Single administrator for both Oracle Database and application software.	Software to control, acquire and process the data with Oracle data base or SQL for high level of data integrity This is deleted


Dr Richa Shri (PI)
Professor

Department of Pharmaceutical Sciences
and Drug Research,
Punjabi University, Patiala